
PRO SL Synus Lift

PRO SL Synus Lift

Ing. C. A. Issoglio & C. S.r.l.
 Via F.lli Di Dio, 68 - 28887 Omegna (VB) - ITALIA
 Tel. +39.0323.828.004 - Fax +39.0323.828.914
 www.resista.it - info@resista.it

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The use of rotating instruments for the Mini-Augmentation Technique through the crest, found its first applications in 2004 with the Dr. F. Cosci Technique.

In the last 18 years this protocol experienced important evolutions both in the technology of materials and engineering of new tools, more effective and less traumatic.

The new generation rotating instruments, with low revs and counterclockwise rotation, minimize the probability of Schneider's membrane perforation without the percussion trauma typical of the osteotomic Summers Technique.

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
The **anticlockwise rotation**, that reverse the involved forces, generates 3 different effects on the implant tunnel perforation, that can revolutionize the known milling logic.

1) Anterior and lateral thrust of the bone, that has been removed by the tip, and of the liquids as well (blood and physiological water).

2) Ejection and "hammer effect" of the drill that produce an improvement in vertical control.

3) Reduction of cutting efficiency in order to protect sensitive anatomical parts.

By pushing the organic material (without removal), an apical bone layer is formed, in front of the drill, that without biomaterial use, is a protection for the lift. The volume of condensed tissue pushes the Schneider membrane with a cushioning effect until it's laterally detached on 360°.

The speed of counterclockwise rotation between to 80 and 100 rpm minimizes the local trauma.

The calibrated stops prevent penetration into the maxillary sinus and allow a gradual millimeter advancement.

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